

PET TUMOR SEGMENTATION: MULTI-OBSERVER VALIDATION OF A GRADIENT-BASED METHOD USING A NSCLC PET PHANTOM



A. D. Nelson¹, M. Werner-Wasik², W. Choi³, Y. Arai⁴, P. F. Faulhaber^{5,6}, N. Ohri², J. W. Piper^{1,7}, K. D. Brockway¹, A. S. Nelson¹

¹MIMvista Corp, Cleveland, OH, ²Thomas Jefferson University Hospital, Philadelphia, PA, ³Beth Israel Medical Center, New York, NY, ⁴University of Pittsburgh Medical Center Health Systems, Pittsburgh, PA, ⁵University Hospitals Case Medical Center, Cleveland, OH, ⁶Case Western Reserve University, Cleveland, OH, ⁷Dept of Computer Science, Wake Forest University, Winston-Salem, NC


Purpose

Consistent and accurate methods for PET tumor segmentation are needed in Radiation Therapy with the growing role of PET for target definition, prognosis, and therapy response assessment. Previously we demonstrated the superior accuracy of a gradient-based PET segmentation method, PET Edge (GRAD), to constant threshold (THRESH) methods using constant activity spherical phantoms¹ which were obtained to emulate clinical activity, acquisition, and reconstruction on five different PET cameras. In a follow-up single observer study we demonstrated the superior accuracy of GRAD compared to constant threshold methods² using published Monte Carlo simulated PET scans of the thorax with "tumors" in commonly encountered locations for lung cancer patients treated in Radiation Therapy³. Our goal in the present study is to determine for multiple observers the accuracy and consistency of GRAD, THRESH and manual contouring (MC) using the same Monte Carlo simulated PET scans of the thorax.

Materials/Methods


Twenty-five realistic digital PET phantoms of the thorax were obtained with 31 simulated tumors of varying size, shape and location. Tumors ranged in size from 7 ml to 264 ml with some tumors surrounded by lung while other tumors were adjacent to the mediastinum or chest wall. Five observers (three Radiation Oncologists, one Radiation Oncologist with Radiology training, and one PET trained Radiologist) segmented each tumor with GRAD, THRESH, and MC. THRESH was performed using thresholds of 25-50% of maximum counts at 5% increments. Tumor volumes for each method were compared to known volumes prior to simulation. Accuracy was measured by calculating the mean absolute % difference for the volume (absVdiff%) and mean % difference for the volume (Vdiff%) for each group using all methods.

Figure 1 Inter-observer Variability




Surrounding normal tissues and heterogeneous tumor activity can result in substantial variability in tumor segmentation. This example demonstrates the robust, consistent results generated using PET Edge compared with threshold and manual contouring.

Figure 2 Phantom Image



For cases where the tumor is in the center of lung (i.e. high source-to-background) and fairly homogeneous, 25% THRESH performs fairly well, however, in the mediastinum with lower source-to-background 25% THRESH performs poorly. GRAD produces more accurate segmentations in both of these scenarios.

Figure 3 Patient Image



Case illustrates limitations of constant threshold in cases of heterogeneous tumor metabolism and areas of decreased source to background ratio (adjacent to chest wall) while GRAD produces more accurate segmentations in both of these scenarios.

Results

Across all observers, GRAD was the most accurate segmentation technique (see table 1) with absVdiff% of 11.3% (12.2% SD). The next most accurate, 25% THRESH, with absVdiff% of 15.1% (16.7%), and MC, with absVdiff% 18.9% (17.3%), were statistically significantly less accurate (p-value < 0.01). GRAD also had the smallest amount of systematic bias (see table 2) with volume differences Vdiff% of -0.45% (16.7%) compared with -7.0% (21.4%) and -15.7% (20.3%) for 25% THRESH and MC, respectively (p-value < 0.01). Finally, GRAD was more consistent between observers than either MC or 25% THRESH in 25/31 and 21/31 lesions, respectively (See Figure 1). Using Levene's test for equality of variances, inter-observer variability was statistically significantly reduced using GRAD compared to MC (p-value < 0.01) but GRAD vs 25% THRESH was non-significant.

Table 1 Volumetric Accuracy

Technique	Observer 1	Observer 2	Observer 3	Observer 4	Observer 5	Overall (Average)
PET Edge	10.979	11.017	13.106	10.783	10.799	11.337
Manual	13.603	25.175	36.541	8.241	10.993	18.911
25%	13.433	11.242	12.776	21.585	16.368	15.081
30%	19.007	17.091	17.198	26.390	21.399	20.217
35%	25.392	23.274	24.389	31.939	25.739	26.147
40%	32.538	31.987	31.849	38.370	34.711	33.891
50%	50.405	49.622	50.332	47.863	51.644	49.973

Mean Absolute Percent Differences in Volume (Vdiff%)

Table 2 Systematic Bias

Technique	Observer 1	Observer 2	Observer 3	Observer 4	Observer 5	Overall (Average)
PET Edge	-0.354	-0.508	-3.870	1.547	0.934	-0.450
Manual	-7.280	-23.373	-35.279	-4.581	-7.964	-15.695
25%	-5.536	-6.389	-4.210	-19.281	0.176	-7.048
30%	-17.156	-12.194	-13.880	-26.248	-12.910	-16.478
35%	-23.680	-21.792	-22.780	-31.939	-24.304	-24.899
40%	-32.538	-29.764	-31.147	-38.370	-34.711	-33.306
50%	-50.405	-49.622	-50.332	-47.862	-51.644	-49.973

Mean Percent Differences in Volume

Conclusions

The gradient-based PET segmentation technique, PET Edge (GRAD) was the most accurate and consistent PET tumor segmentation technique across observers. GRAD has the potential to play an important role in Radiation Therapy through more accurate and consistent PET Tumor segmentation.

References

- Shen G, Nelson D, Adler L. PET Tumor Segmentation: Comparison of Gradient-Based Algorithm to Constant Threshold Algorithm. Med Phys. 2007;34:2395.
- Nelson D, et al. PET Tumor Segmentation: validation of a gradient-based method using a NSCLC PET Phantom. J Nucl Med. 2009; 50 (Supplement 2):1659.
- Aristophanous M, Penney BC, Pelizzari CA. The development and testing of a digital PET phantom for the evaluation of tumor volume segmentation techniques. Med Phys. 2008;35:3331-3342.