



**mim**  
**VIEWER™**

Comparison Chart

MIMviewer™

MIMviewer™  
Pro

**MIMviewer Features List**

DICOM Viewing on a Remote Computer	X	X
Display all PET/CT and most SPECT/CT, PET, CT, MR, SPECT	X	X
Display of SUV and Hounsfield Units (Voxel)	X	X
Incremental Volume Blending	X	X
Rotating MIP Movie with/without Depth Shading	X	X
Display Pre-aligned Fused Images	X	X
Patient Information	X	X
Multi-Planar Reconstruction (MPR) Display	X	X
Customizable Display of Single and Fused Volumes	X	X
Image Display Tools	X	X
Slice Selection	X	X
Triangulation	X	X
Magnification	X	X
Contrast	X	X
Color Tables	X	X
CT Preset Tissue Densities	X	X
Screen Capture and Print Support	X	X
Structured DICOM Report	X	X
Patient Anonymization	X	X
User Configuration of General and Interface Parameters	X	X
Multiple Language Support ( <i>English, French, Spanish, Chinese, Japanese</i> )	X	X
Windows XP, Vista, and Mac OSX 10.4+	X	X
Full Volume Blending		X
3D PET SUV Tool		X
Zoom & Pan		X
Linear Measure Tool		X
DICOM Store/Receive & Query/Retrieve		X
Generate MIMviewer DICOM CD		X